

BioElectronics

bio-impedance measurement and modeling

F. Kölbl

Maître de Conférences

florian.kolbl@enseirb-matmeca.fr



Why modeling?



Scientific aphorism

... all models are approximations. Essentially,
all models are wrong, but some are useful. However, the
approximate nature of the model must always be borne in mind...
George Box, *Empirical Model-Building and Response Surfaces*, 1987

In a largely pluri-disciplinary, context a model is an abstraction that we can discuss, whatever our field of expertise. It is an opportunity for engineers to meet biology

let's discuss bioelectronics with it!

Today's agenda


Lecture

- What are the basic physics of (passive) bioelectricity?
- How to measure and model bio-impedance?

Lab


- Electrodes and tissue impedance measurement,
- Bioelectrical measurement and modeling,
- Bioelectronics and potatoes!

Bioelectric interfaces



Malmivuo, J., Plonsey, R. (1995). *Bioelectromagnetism: principles and applications of bioelectric and biomagnetic fields*. Oxford University Press, USA.

The cell






Living cell

micrometric machine block of a living organism, with:

- chemical, molecular and protidic capabilities,
- procedure storage capabilities,
- potential electrical activity, at least electrical properties.

Variety of cells



muscle cells:

- mechanically active,
- electro sensitive,

neurons:

- electrically active,
- electro sensitive,

bone cells:

- not excitable,
- have passive properties

All cells have common passive electrical and dielectrical properties, tissue-impedance is a singular characteristic.

1 Tissue passive properties

- Conduction in ionic media
- Unveiling the membrane

2 Bio-impedance measurement

- Connecting to the tissues: electrodes
- Properties measurement: how-to
- Example application

Electrical Conduction in biological tissues

Electrical conduction has a different nature considering the medium

Electrical Circuits

charge carrier: electron

elec. charge:

$$-1 \cdot e$$

Tissue

charge carrier: ions

elec. charges:

$$\begin{aligned} Na^+, K^+ &: +1 \cdot e \\ Cl^-, HCO_3^- &: -1 \cdot e \\ Ca^{2+}, Mg^{2+} &: +2 \cdot e \dots \end{aligned}$$

current:

$$i = \frac{dQ}{dt}$$


current:

$$i = \sum_{ions} I_{ion}$$

each ion can move due to *migration, diffusion, convection*

Resistance, conductance

In a homogeneous material:



electrical resistance given by:

$$R = \rho \frac{l}{S} = \frac{1}{\sigma} \cdot \frac{l}{S}$$

with ρ the resistivity in $\Omega \cdot m$ or σ the conductivity in $S \cdot m^{-1}$,

material	conductivity ($S \cdot m^{-1}$)
coper	$6 \cdot 10^7$
germanium	2.17
deionized water	$5.5 \cdot 10^{-6}$

Resistivity, conductivity

In a ionic solution:

the conductivity is given by:

$$\sigma = \sum_{k \text{ ions}} \Lambda_k c_k$$

with c the chemical concentration in $mol \cdot L^{-1}$ and Λ the molar conductivity in $S \cdot m^2 \cdot mol^{-1}$

and in strong electrolytes, at very low concentration, as in living organisms, $\Lambda \approx \Lambda_0$ independant from the concentration.

Cation	Λ_0 in $S \cdot cm^2 \cdot mol^{-1}$	Anion	Λ_0 in $S \cdot cm^2 \cdot mol^{-1}$
H^+ / H_3O^+	350	OH^-	198
Na^+	50	Cl^-	76
K^+	74	HCO_3^-	45
Ca^{2+}	119	CO_3^-	72

Example: Conductivity of the 0.9% Saline solution

9g of NaCl per Liter of water
the atomic mass of NaCl is $58.5 \text{ g} \cdot \text{mol}^{-1}$



recall $\Lambda_{0,\text{Na}^+} = 50 \text{ S} \cdot \text{cm}^2 \cdot \text{mol}^{-1}$ and $\Lambda_{0,\text{Cl}^-} = 76 \text{ S} \cdot \text{cm}^2 \cdot \text{mol}^{-1}$

$$[\text{Na}^+] = [\text{Cl}^-] = \frac{m_{\text{NaCl}}}{M_{\text{NaCl}}} = \frac{9}{58.5} = 0.154 \text{ mol} \cdot L^{-1}$$



$$\sigma = \Lambda_{0,\text{Na}^+} [\text{Na}^+] + \Lambda_{0,\text{Cl}^-} [\text{Cl}^-] = \frac{0.154 (50 + 76)}{1000} \approx 19 \text{ mS} \cdot \text{cm}^{-1}$$


(divided by 1000 to convert the L^{-1} in cm^{-3})

Cell membranes

Physical separation between the intra- and extra-cellular medium

- about 5nm thick
- phospholipid-bilayer: one layer is composed of one hydrophobic and one hydrophilic lipid that self assemble in membrane





*Lipidic (insulating) membrane, separating two conductive electrolytes
that ionic moving charges cannot cross
⇒ Equivalent to a capacitance*

Computation of the membrane capacitance

Data

for a 5 nm thick membrane,
 $\epsilon_0 = 8.85418782 \cdot 10^{-12} \text{ m}^{-3} \cdot \text{kg}^{-1} \cdot \text{s}^4 \cdot \text{A}^2$
 the relative membrane permittivity is $\epsilon_r = 5$

$$C = \epsilon_0 \epsilon_r \frac{S}{e} = \tilde{c}_M \cdot S$$

where \tilde{c}_M is the specific membrane capacity ($\text{F} \cdot \text{m}^{-2}$)

$$\tilde{c}_M = \frac{\epsilon_0 \epsilon_r}{e} = \frac{5 \times 8.85418782 \cdot 10^{-12}}{5 \cdot 10^{-9}} \approx 8.85 \text{ mF} \cdot \text{m}^{-2} = 0.885 \mu\text{F} \cdot \text{cm}^{-2}$$

$\tilde{c}_M = 1 \mu\text{F} \cdot \text{cm}^{-2}$ is a common value in the literature

Suming up passive properties


Intra/extracellular medium

- two resistive media
- small ionic concentration changes enable to consider it as constant resistivity, (especially in extra-cellular space)

All cellular membranes

capacitive

first approximation tissue impedance model:



warning: tissue only, no electrode

outline


1 Tissue passive properties

- Conduction in ionic media
- Unveiling the membrane


2 Bio-impedance measurement

- Connecting to the tissues: electrodes
- Properties measurement: how-to
- Example application

Electrode classification



First schematic view



Not a direct electrical acces to the tissue
One or more materials directly in contact with the tissue

Which material? (1/2)

		toxicity	reactivity
conductors	Gold	non-toxic	non-reactive
	Silver	toxic	
	Copper	toxic	
	Iron	toxic	
	Stainless Steel	non-toxic	
	Platinum	toxic	
	Tantalum		reactive
	Titanium		
	Tungsten		non-reactive
	Gold-nickel-chromium	non-toxic	
	Gold-palladium-rhodium	non-toxic	
	Nickel-chromium (Nichrome)	non-toxic	reactive
	Nickel-chromium-molybdenum	non-toxic	
	Nickel-titanium (Nitinol)		
	Platinum-iridium	non-toxic	
	Platinum-nickel	non-toxic	
	Platinum-rhodium	non-toxic	
	Platinum-tungsten	non-toxic	
	Platinized platinum (Pt black)	non-toxic	
	...		

Which material? (2/2)


		toxicity	reactivity
Semi-conductors	...		
	Silicon		
	Germaniom	toxic	non-reactive
			biocompatible
Insulators	Alumina ceramic		non-reactive
	Araldite (epoxy plastic resin)		reactive
	Polyethylene		non-reactive
	Polyimide		
	Polypropylene		non-reactive
	Silicon dioxide (Pyrex)		reactive
	Teflon TFE (high purity)		non-reactive
	Teflon TFE (shrinkable)		reactive
	Titanium dioxide		reactive

adapted from Merrill, D. R., Bikson, M., Jefferys, J. G. (2005). Electrical stimulation of excitable tissue: design of efficacious and safe protocols. Journal of neuroscience methods, 141(2), 171-198.

commonly used materials in electronics and micro-electronics (gold, stainless steel, silicon, polyimide) can be used!

warning: no copper

The double layer (1/2)



The double layer (2/2)

Helmholtz Capacitance

$$C_H = \epsilon_0 \epsilon_r \frac{A}{d_H}$$

- d_H is a constant
- A the effective electrode Surface Area
about $230 \mu F \cdot cm^{-2}$
value depends on surface roughness

Gouy-Chapman Capacitance

$$C_D = \frac{\epsilon_0 \epsilon_r}{L_D} \cosh \left(\frac{q_i \Phi_0}{2RT} \right)$$

- q_i is the ion charge
- $L_D = \sqrt{\frac{\epsilon_0 \epsilon_r}{2RT c_i q_i}}$, with c_i the ion concentration
- with Φ_0 the junction voltage
about $50 \mu F \cdot cm^{-2}$
possibly (voltage) non-linear

overall specific capacitance value about $40 \mu F \cdot cm^{-2}$

At the junction between a metal and a conductive electrolyte: electrical voltage
(Electrochemical half-cell potential) depending on the metal

Material	Reaction	Potential
Aluminium	$Al^{3+} + 3e^-$	-1.67V
Iron	$Fe^{2+} + 2e^-$	-0.441V
Silver	$Ag^+ + e^-$	+1.7996V
Platinum	$Pt^{2+} + 2e^-$	+1.2V
Gold	$Au^{3+} + 3e^-$	+1.52V
	$Au^+ + e^-$	+1.83V
H_2	$2H^+ + 2e^-$	0.000V (Reference)

at $T = 298K$

Note that if symmetrical materials \rightarrow overall voltage = 0V


The Faradaic impedance

charges may be shared by redox reactions:
transfer of electrons between the two phases (metal, electrolyte)

- Reactions depends on the material,
- highly (voltage) non-linear,
- complex modeling (*resistive but not that much, nor capacitive...*),
- in electrochemistry, considered as a **Constant Phase Element**.



let us call it Z_F , we will speak about it later

First approximation physical model







Impedance Measurement: 2-points configuration:

For bio-potential recording and impedance sensing
Electrical Impedance Spectroscopy







2-points: electrode size




2-points: parasitics

- Larg
- Low
- Sma
- Higl



4-points configuration:

For bio-potential impedance sensing only




F. Kölbl (ENSEIRB-Matmeca)

BioElectronics

29 / 37

Physiological reactions




Proteic layer (seconds) → Cell migration (minutes) → Growth (hours) → Differentiation (days) → Matrix secretion (weeks)

Brette, R., Destexhe, A. (2012). Handbook of neural activity measurement. Cambridge University Press.

Problems after implantation directly on electrodes


pacemaker implantation at time $t = 0$



explantation after $t = 5$ years




Mechanism of fibrosis



WICK, Georg, GRUNDTMAN, Cecilia, MAYERL, Christina, et al. The immunology of fibrosis. Annual review of immunology, 2013, vol. 31, p. 107-135.


A change in the medium



Warning: the reaction is highly local and unpredictable

Polikov, V. S., Tresco, P. A., Reichert, W. M. (2005). Response of brain tissue to chronically implanted neural electrodes. *Journal of neuroscience methods*, 148(1), 1-18.

Electrical consequences





Williams, J. C., Hippensteel, J. A., Dilgen, J., Shain, W., Kipke, D. R. (2007). Complex impedance spectroscopy for monitoring tissue responses to inserted neural implants. *Journal of neural engineering*, 4(4), 410.


However: no clear correlation between impedance change and physiological reaction

Example of investigation (ETIS work)


Investigation to correlate physiology (*biology*) and impedance (*electronics*)



design of an *in vitro* measurement bench on gold-microelectrodes





impedance control



tissu evolution evaluation
Electrode
Protein fibers
Cell Nucleus

Work of Edwin De Roux - ETIS - 2018


Cell proliferation monitoring



Correlation between proliferation, apoptosis and impedance
biomarker for large populations monitoring


impedance model :

Ex vivo experiments



work of A. Degache, N. Lewis (IMS), O. Bernus (IHU Lyric), F. Kolbl (ETIS)

$$Z(j\omega) = R \frac{1 + \left(\frac{jf}{f_\alpha}\right)^\alpha}{\left(\frac{jf}{f_\alpha}\right)^\alpha} \frac{\left(\frac{jf}{f_\beta}\right)^\beta}{1 + \left(\frac{jf}{f_\beta}\right)^\beta}$$



first demonstration of the possibility to use electrode as a sensor to discriminate tissues